
Analytical Fittings

Union Connector – Low Dead Volume

For fractional tube





CPI™ PART NO.	A-LOK® PART NO.	INTER-CHANGES WITH	INCHES													INTERNAL VOLUME
			T_1 TUBE O.D.	T_2 TUBE O.D.	†A	†B	†C	D_1	D_2	E INTERNAL OPENING	G_1	G_2	H	W HEX		
1-1 Z7HBZ7-SS	1-1 Z7HLZ7	IFO-6GC	1/16	1/16	1.26	.21	.21	.41	.41	.013	.25	.25	.84	1/4		8.7 x 10-5cc
2-1 Z7HBZ7-SS	2-1 Z7HLZ7	–	1/8	1/16	1.53	.31	.21	.56	.41	.013	.38	.25	1.02	7/16		8.7 x 10-5cc
2-2 Z7HBZ7-SS	2-2 Z7HLZ7	–	1/8	1/8	1.81	.31	.31	.56	.56	.052	.38	.38	1.19	7/16		9.7 x 10-2cc

†Average Value

Dimensions for reference only, subject to change.